



**ELECTRO
SURGICAL**



CATALOGUE
Supplier of Quality Electro Surgical Instruments

www.propasintl.com.pk

PRO-PAS INTERNATIONAL COMPANY PROFILE

Pro-Pas International Electro Surgical is a leading innovator in the field of electrosurgical equipment and medical devices, providing state-of-the-art solutions for surgical procedures across multiple medical disciplines. Our company is dedicated to enhancing surgical precision, patient safety, and operational efficiency through advanced technology and innovation.

With a strong commitment to quality and excellence, Electro Surgical has earned a reputation for delivering high-performance electrosurgical generators, cutting-edge instruments, and reliable accessories that meet the needs of modern healthcare institutions. Our products are designed with ergonomic efficiency, safety features, and seamless integration, ensuring optimal results in both routine and complex surgical procedures.

OUR MISSION

To empower medical professionals with cutting-edge electrosurgical technology that enhances surgical efficiency, reduces risks, and improves patient outcomes.

OUR VISION

To be a global leader in electrosurgical innovation, setting new standards in medical technology and ensuring accessibility to superior surgical equipment.

TABLE OF CONTENTS


1	<i>European Bipolar Forceps</i>	1 to 8
2	<i>European Irrigation Bipolar Forceps</i>	9
3	<i>European Non Stick Bipolar Forceps</i>	10 to 13
4	<i>Silver Clide Non Stick Bipolar Forceps</i>	14 to 20
5	<i>Us Bipolar Forceps</i>	21 to 25
6	<i>Monopolar Forceps</i>	26 to 27
7	<i>Us Irrigation Bipolar Forceps</i>	28 to 29
8	<i>Us Non Stick Bipolar Forceps</i>	30 to 33
9	<i>Hand Switch Monopolar Forceps</i>	34
10	<i>Diathermy Instruments</i>	35 to 37
11	<i>Disposable Bipolar Forceps</i>	38 to 43
12	<i>Ob/Gyn Leep Instruments</i>	44 to 47
13	<i>Laparoscopic & Arthoscopic Electrode</i>	48 to 49
14	<i>Electrosurgical Pancil & Electrode</i>	50 to 55
15	<i>Electrosurgical Cable</i>	56 to 59
16	<i>Bipolar Artery Sealer</i>	60
17	<i>Bipolar Scissors</i>	61

The Pro-Pas electro bipolar forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable Bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.

Advantages of Reusable European Bipolar Forceps

- 1- Highly polished tips reduce tissue adherence.*
- 2- Light closing force causes less fatigue over time.*
- 3- Nylon coating resists abrasion and provides superb insulation.*
- 4- Reusable, durable and designed for multiple uses.*
- 5- Precision tip alignment aids fine control.*
- 6- The Instruments is autoclavable.*
- 7- Compatible with all common electrosurgery unit.*






100




101



102




103



104




105




Ref	Item Name	Tip Size (mm)	Length Size (cm)
100	MCPHERSON	0.5	10.2
101	MCPHERSON ANGLED	0.5	10.2
102	JEWELLER #5	0.5, 0.7, 1.0	11.5
103	JEWELLER #5 ANGLED	0.5, 0.7, 1.0	11.5
104	JEWELLER #7	0.5, 0.7, 1.0	11.5
105	ADSON	0.5, 0.7, 1.0	12




Ref	Item Name	Tip Size (mm)	Length Size (cm)
106	ADSON	0.5, 0.7, 1.0	15
107	STRAIGHT	0.5, 1.0, 2.0	15
108	STRAIGHT	0.5, 1.0, 2.0	18
109	STRAIGHT	0.5, 1.0, 2.0	20
110	STRAIGHT	0.5, 1.0, 2.0	22
111	STRAIGHT	0.5, 1.0, 2.0	24







Ref	Item Name	Tip Size (mm)	Length Size (cm)
112	STRAIGHT ANGLED	0.5, 1.0, 2.0	24
113	STRAIGHT ANGLED	0.5, 1.0, 2.0	22
114	STRAIGHT ANGLED	0.5, 1.0, 2.0	20
115	STRAIGHT ANGLED	0.5, 1.0, 2.0	18
116	STRAIGHT ANGLED	0.5, 1.0, 2.0	15




Ref	Item Name	Tip Size (mm)	Length Size (cm)
117	BAYONET	0.5, 1.0, 2.0	24
118	BAYONET	0.5, 1.0, 2.0	22
119	BAYONET	0.5, 1.0, 2.0	20
120	BAYONET	0.5, 1.0, 2.0	18
121	BAYONET	0.5, 1.0, 2.0	16.5






Ref	Item Name	Tip Size (mm)	Length Size (cm)
122	BAYONET ANGLED UP	0.5, 1.0, 2.0	24
123	BAYONET ANGLED UP	0.5, 1.0, 2.0	22
124	BAYONET ANGLED UP	0.5, 1.0, 2.0	20
125	BAYONET ANGLED UP	0.5, 1.0, 2.0	18
126	BAYONET ANGLED UP	0.5, 1.0, 2.0	16.5



Ref	Item Name	Tip Size (mm)	Length Size (cm)
128	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	24
130	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	22
132	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	20
134	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	18
136	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	16.5




It is well known that depending on the generator and the surface quality of the forceps jaws, coagulated tissue tends to more or less stick to the instruments jaw. This unwanted effect can occur particularly oftne with liver tissue and is critical in neurosurgery. For these applications in particular, we have designed our Bipolar Forceps with integrated irrigation channedl. The coagulation area is irrigated through a built-in irrigation channedl inside the forceps branch which can be connected to a standard irrigation hose. By means of this irrigation the sticking of the tissue is substantially reduced and surgical work is made easier.


Advantages of Reusable European Irrigation Bipolar Forceps

- 1-** *A universal hub connector allows irrigation tubing to be connected easily.*
- 2-** *The ergonomically designed handle provi es less fatigue and a well balanced instruments.*
- 3-** *Nylon coating resists abrasion and provides superb insulation.*
- 4-** *For high visibility our irrigation tube can be recessed into the forceps blank (buried irrigation).*
- 5-** *Compatible with all common electrosurgery unit.*






Ref	Item Name	Tip Size (mm)	Length Size (cm)
138	STRAIGHT IRRIGATION	0.5, 1.0, 2.0	24
140	STRAIGHT IRRIGATION	0.5, 1.0, 2.0	22
142	STRAIGHT IRRIGATION	0.5, 1.0, 2.0	20
144	STRAIGHT IRRIGATION	0.5, 1.0, 2.0	18
146	STRAIGHT IRRIGATION	0.5, 1.0, 2.0	15





Tip Chart




Ref	Item Name	Tip Size (mm)	Length Size (cm)
148	BAYONET IRRIGATION	0.5, 1.0, 2.0	24
150	BAYONET IRRIGATION	0.5, 1.0, 2.0	22
152	BAYONET IRRIGATION	0.5, 1.0, 2.0	20
154	BAYONET IRRIGATION	0.5, 1.0, 2.0	18
156	BAYONET IRRIGATION	0.5, 1.0, 2.0	16.5

During the caogulation process, tissue and tips of forceps can stick together. Pro-Pas non stick bipolar forceps are made with a special silver alloy and therefore prevent tissue from sticking. Disturbing and time consuming disruptions for cleaning the tips during surgery are thus avoided. Due to the good thermal qualities of silver the termperature of the forceps tips does not exceed 100'C during coagulation, which enhances the non-stick effect. With over ten years of experience we can guarantee that out bipolar forceps do not stick. The non-stick effect remains constant and is not reduced by frequent sterilization. Our non-stick bipolar forceps together with the appropriate connection cables are compatible with all common types of high-frequency units such as Erbe, Martin, Berchtold, Valleylab. Pro-Pas Bipolar forceps non-stick are available in straight, angled or bayonet shape. Available lengths range from 110 to 300 mm. Tips can be made from 0.25 to 2.0 mm width. Forceps are also available with different measureme' ts as well as with flat pin or 2-pin plug.

Advantages of Reusable Non-Stick European Bipolar Forceps

- 1-** Ten years of experience speak for themselves. Our forceps do not stick-never!
- 2-** No cleaning during surgery.
- 3-** Especially high thermal qualities.
- 4-** All Non-Stick Bipolar Forceps are re-usable and autoclavable.
- 5-** Compatible with all common electrosurgerv unit.





Ref	Item Name	Tip Size (mm)	Length Size (cm)
158	MCPHERSON	0.5, 1.0,	10.2
160	MCPHERSON ANGLED	0.5, 1.0,	10.2
162	JEWELLER #5	0.5, 1.0,	11.5
164	JEWELLER #5 ANGLED	0.5, 1.0,	11.5
166	ADSON	0.5, 1.0,	12



168



170



172




174



176


Ref	Item Name	Tip Size (mm)	Length Size (cm)
168	STRAIGHT	0.5, 1.0, 2.0	24
170	STRAIGHT	0.5, 1.0, 2.0	22
172	STRAIGHT	0.5, 1.0, 2.0	20
174	STRAIGHT	0.5, 1.0, 2.0	18
176	STRAIGHT	0.5, 1.0, 2.0	15

Tip Chart





Tip Chart




Ref	Item Name	Tip Size (mm)	Length Size (cm)
178	BAYONET	0.5, 1.0, 2.0	24
180	BAYONET	0.5, 1.0, 2.0	22
182	BAYONET	0.5, 1.0, 2.0	20
184	BAYONET	0.5, 1.0, 2.0	18
186	BAYONET	0.5, 1.0, 2.0	16.5

The Pro-Pas bipolar forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable Bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.

Advantages of Reusable US Bipolar Forceps

- 1- Highly polished tips reduce tissue adherence.*
- 2- Light closing force causes less fatigue over time.*
- 3- Nylon coating resists abrasion and provides superb insulation.*
- 4- Reusable, durable and designed for multiple uses.*
- 5- Precision tip alignment aids fine control.*
- 6- The Instruments is autoclavable.*
- 7- Compatible with all common electrosurgery unit.*







Ref	Item Name	Tip Size (mm)	Length Size (cm)
188	JEWELLER #5	0.5, 1.0,	11.5
190	ADSON	0.5, 1.0,	12

SILVER GLIDE NON STICK BIPOLAR FORCEPS



Ref	Item Name	Tip Size (mm)	Length Size (cm)
192	STRAIGHT	0.5, 1.0, 2.0	24
194	STRAIGHT	0.5, 1.0, 2.0	22
196	STRAIGHT	0.5, 1.0, 2.0	20
198	STRAIGHT	0.5, 1.0, 2.0	18
200	STRAIGHT	0.5, 1.0, 2.0	15







Ref	Item Name	Tip Size (mm)	Length Size (cm)
202	STRAIGHT ANGLED	0.5, 1.0, 2.0	24
204	STRAIGHT ANGLED	0.5, 1.0, 2.0	22
206	STRAIGHT ANGLED	0.5, 1.0, 2.0	20
208	STRAIGHT ANGLED	0.5, 1.0, 2.0	18
210	STRAIGHT ANGLED	0.5, 1.0, 2.0	15




Ref	Item Name	Tip Size (mm)	Length Size (cm)
212	BAYONET	0.5, 1.0, 2.0	24
214	BAYONET	0.5, 1.0, 2.0	22
216	BAYONET	0.5, 1.0, 2.0	20
218	BAYONET	0.5, 1.0, 2.0	18
220	BAYONET	0.5, 1.0, 2.0	16.5







Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
222	BAYONET ANGLED UP	0.5, 1.0, 2.0	24
224	BAYONET ANGLED UP	0.5, 1.0, 2.0	22
226	BAYONET ANGLED UP	0.5, 1.0, 2.0	20
228	BAYONET ANGLED UP	0.5, 1.0, 2.0	18
230	BAYONET ANGLED UP	0.5, 1.0, 2.0	16.5



Ref	Item Name	Tip Size (mm)	Length Size (cm)
232	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	24
234	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	22
236	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	20
238	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	18
240	BAYONET ANGLED DOWN	0.5, 1.0, 2.0	16.5




The Pro-Pas bipolar forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable Bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.


Advantages of Reusable US Bipolar Forceps

- 1- Highly polished tips reduce tissue adherence.*
- 2- Light closing force causes less fatigue over time.*
- 3- Nylon coating resists abrasion and provides superb insulation.*
- 4- Reusable, durable and designed for multiple uses.*
- 5- Precision tip alignment aids fine control.*
- 6- The Instruments is autoclavable.*
- 7- Compatible with all common electrosurgery unit.*






Ref	Item Name	Tip Size (mm)	Length Size (cm)
242	MCPHERSON	0.5, 0.7, 1.0	10.2
244	MCPHERSON ANGLED	0.5, 0.7, 1.0	10.2
246	JEWELLER #5	0.5, 0.7, 1.0	11.5
248	JEWELLER #4	0.5, 0.7, 1.0	11.5
250	IRIS STRAIGHT	0.5, 0.7, 1.0	10.8






Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
252	IRIS ANGLED	0.5, 0.7, 1.0	10.8
254	ADSON	0.5, 0.7, 1.0	12
256	SEMKIN	0.5, 0.7, 1.0	14
258	CUSHING	0.5, 0.7, 1.0, 1.5	15.9
260	CUSHING	0.5, 0.7, 1.0, 1.5	17.8
262	GERALD STRAIGHT	0.5, 1.0, 1.5	17.8
264	GERALD BAYONET	0.5, 1.0, 1.5	18.4




Ref	Item Name	Tip Size (mm)	Length Size (cm)
266	JANSEN BAYONET	0.5, 1.0, 1.5, 2.0	20.4
268	JANSEN BAYONET	0.5, 1.0, 1.5	17.8
270	JANSEN BAYONET	0.5, 1.0, 1.5	15.2





Tip Chart





Ref	Item Name	Tip Size (mm)	Length Size (cm)
272	MALIS	0.5, 1.0, 1.5	17.8
274	SCOVILLE-GREENWOOD	0.5, 1.0, 1.5	19.1
276	YASARGIL	0.5, 1.0, 1.5	19.7
278	HARDY BAYONET	0.5, 1.0, 1.5	21.6

Our Monopolar instruments are tested at high volts to ensure proper functioning during the most stringent operations. With the best quality stainless steel, excellent finishing and superior insulation, we provide worldwide guarantee for all our Monopolar instruments.

Advantages of Reusable Monopolar Forceps

- 1- Lightweight and well balanced to prevent fatigue.*
- 2- Highly polished tips minimise sticking.*
- 3- Nylon coating resists abrasion and provides superb insulation.*
- 4- Precision tip alignment aids fine control.*
- 5- Compatible with all common electrosurgery unit.*






Ref	Item Name	Tip Size (mm)	Length Size (cm)
280	Jeweller #4	Needle	11.5
282	Jeweller #7	Needle	11.5
284	Adson	1.0	12.0
286	Potts-Smith	2.0	17.8
288	Potts-Smith	2.0	21.0
290	Potts-Smith	2.0	24.1
292	Jansen Bayonet	1.5	13.3
294	Jansen Bayonet	1.5	15.9
296	Jansen Bayonet	2.0	19.1
298	Jansen Bayonet	2.0	21.0
300	Jansen Bayonet UP	2.0	19.1
302	Jansen Bayonet UP	2.0	21.0

It is well known that depending on the generator and the surface quality of the forceps jaws, coagulated tissue tends to more or less stick to the instruments jaw. This unwanted effect can occur particularly often with liver tissue and is critical in neurosurgery. For these applications in particular, we have designed our Bipolar Forceps with integrated irrigation channel. The coagulation area is irrigated through a built-in irrigation channel inside the forceps branch which can be connected to a standard irrigation hose. By means of this irrigation the sticking of the tissue is substantially reduced and surgical work is made easier.

Advantages of Reusable US Irrigation Forceps

- 1-** *A universal hub connector allows irrigation tubing to be connected easily.*
- 2-** *The ergonomically designed handle provides less fatigue and a well balanced instruments.*
- 3-** *For high visibility our irrigation tube can be recessed into the forceps blank (buried irrigation).*
- 4-** *Compatible with all common electrosurgery unit.*





Ref	Item Name	Tip Size (mm)	Length Size (cm)
304	STRAIGHT	0.5, 1.0, 1.5	15.9
306	GERALD BAYONET	0.5, 1.0, 1.5	18.4
308	JANSEN BAYONET	0.5, 1.0, 1.5	19.1
310	HARDY BAYONET	0.5, 1.0, 1.5	21.6

During the caogulation process, tissue and tips of forceps can stick together. Pro-Pas non stick bipolar forceps are made with a special silver alloy and therefore prevent tissue from sticking. Disturbing and tims consuming disruptions for cleaning the tips during surgery are thus avoided. Due to the good thermal qualities of silver the termperature of the forceps tips does not exceed 100C during coagulation, which enhances the non stick effect. With over ten years of experience we can guarantee that our bipolar forceps do not stick. The non stick effect remains constant and is not reduced by frequent sterilization. Our non stick bipolar forceps together with the appropriate connection cables are compatible with all common types of high-frequency units such as Erbe, Martin, Berchtold, Valleylab. Pro-Pas Bipolar forceps non stick are available in straight, angled or bayonet shape. Available lengths range from 110 to 300 mm. Tips can be made from 0.25 to 2.0 mm width. Forceps are also available with different measurements as well as with flat pin or 2-pin plug.

Advantages of Rusable Non-Stick European Bipolar Forceps

- 1-** *Ten years of experience speak for themselves. Our forceps do not stick-never!*
- 2-** *No cleaning during surgery.*
- 3-** *Especially high thermal qualities.*
- 4-** *All Non-stick Bipolar Forceps are re-usable and autooclavable.*
- 5-** *Compatible with all common electrosurgery unit.*





312



314



316



318




320



322


Tip Chart





Ref	Item Name	Tip Size (mm)	Length Size (cm)
312	MCPHERSON	NEEDLE, 0.5	10.2
314	MCPHERSON ANGLED	NEEDLE, 0.5	10.2
316	JEWELLER# 5 ANGLED	NEEDLE, 0.5	11.5
318	JEWELLER# 4	NEEDLE, 0.5	11.5
320	ADSON	0.5, 1.0	12
322	SEMKIN	0.5, 1.0	14



Ref	Item Name	Tip Size (mm)	Length Size (cm)
324	CUSHING	0.5, 1.0, 1.5	17.8
326	CUSHING	0.5, 1.0, 1.5	15.9
328	GERALD BAYONET	0.5, 1.0, 1.5	18.4
330	JANSEN BAYONET	0.5, 1.0, 1.5	18
332	ADSON	1.0, 1.5	19.7






Ref	Item Name	Tip Size (mm)	Length Size (cm)
334	Jeweller #4	0.5	11.5
336	Semkin	0.7	14.0
338	Jansen Bayonet	0.7	19.1
340	Scoville-Greenwood	1.5	19.7



Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
342	Jansen Bayonet	2.0	19.0
344	Jansen Bayonet	2.0	21.0
346	Potts-Smith	2.0	21.0
348	Potts-Smith	2.0	24.0

Our Diathermy Instruments are tested at high volts to ensure proper functioning during the most stringent operations. With the best quality stainless steel, excellent finishing and superior insulation, we provide worldwide guarantee for all our Diathermy Instruments. Diathermy Instruments are available in four base types. Models popular Diathermy Instruments: Adson, Jeweller, McIndoe, Standard, Bayonet, Artery Forceps, Wilson Hey Forceps, Diathermy Scissors etc in various lengths, tip sizes and 4 plug types.

Advantages of Reusable Diathermy Instruments

- 1-** *Lightweight and well balanced to prevent fatigue.*
- 2-** *Highly polished tips minimise sticking.*
- 3-** *Precision tip alignment aids fine control.*
- 4-** *Compatible with all common electrosergery unit.*





ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
350	Tuner Warwick Scissors (STR)	-	17.8
352	Tuner Warwick Scissors (CRV)	-	17.8
354	Tuner Warwick Scissors (STR)	-	23.0
356	Tuner Warwick Scissors (CRV)	-	23.0
358	Dressing Scissors (STR)	-	17.8
360	Dressing Scissors (CRV)	-	17.8
362	Dressing Scissors (STR)	-	23.0
364	Dressing Scissors (CRV)	-	23.0
366	Wilson Hay Forceps (STR)	-	17.8
368	Wilson Hay Forceps (CRV)	-	17.8
370	Wilson Hay Forceps (STR)	-	23.0
372	Wilson Hay Forceps (CRV)	-	23.0
374	Artery Forceps (STR)	-	17.8
376	Artery Forceps (CRV)	-	17.8
378	Artery Forceps (STR)	-	23.0
380	Artery Forceps (CRV)	-	23.0
382	Abbey Needle	-	17.5
384	Richie Forceps	-	17.8
386	Richie Forceps	-	20.0
388	Richie Forceps	-	23.0

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
390	Tuner Warwick Scissors (STR)	-	17.8
392	Tuner Warwick Scissors (CRV)	-	17.8
394	Tuner Warwick Scissors (STR)	-	23.0
396	Tuner Warwick Scissors (CRV)	-	23.0
398	Dressing Scissors (STR)	-	17.8
400	Dressing Scissors (CRV)	-	17.8
402	Dressing Scissors (STR)	-	23.0
404	Dressing Scissors (CRV)	-	23.0
406	Wilson Hay Forceps (STR)	-	17.8
408	Wilson Hay Forceps (CRV)	-	17.8
410	Wilson Hay Forceps (STR)	-	23.0
412	Wilson Hay Forceps (CRV)	-	23.0
414	Artery Forceps (STR)	-	17.8
416	Artery Forceps (CRV)	-	17.8
418	Artery Forceps (STR)	-	23.0
420	Artery Forceps (CRV)	-	23.0
422	Abbey Needle	-	17.5
424	Richie Forceps	-	17.8
426	Richie Forceps	-	20.0
428	Richie Forceps	-	23.0



Adson

430



McINDOE

432



Gillies

434



Dissecting

442



Gerald

450



Adson Bayonet

452



Adson Bayonet UP

456



Remington Hobb

462

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
430	Adson	1.0	12.0
432	McINDOE	1.0, 2.0	15.0
434	Gillies	-	15.0
436	Gillies	-	17.8
438	Gillies	-	20.0
440	Gillies	-	24.0
442	Dissecting	1.0, 2.0	15.0
444	Dissecting	1.0, 2.0	17.8
446	Dissecting	1.0, 2.0	21.0
448	Dissecting	1.0, 2.0	24.1
450	Gerald 1x2	-	17.8
452	Adson Bayonet	1.0, 2.0	16.0
454	Adson Bayonet	1.0, 2.0	19.1
456	Adson Bayonet Angled UP	1.0, 2.0	16.0
458	Adson Bayonet Angled UP	1.0, 2.0	19.1
460	Adson Bayonet Angled UP	1.0, 2.0	21.0
462	Remington Hobb	1.0, 2.0	24.1

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
464	Adson	1.0	12.0
466	McINDOE	1.0, 2.0	15.0
468	Gillies	-	15.0
470	Gillies	-	17.8
472	Gillies	-	20.0
474	Gillies	-	24.0
476	Dissecting	1.0, 2.0	15.0
478	Dissecting	1.0, 2.0	17.8
480	Dissecting	1.0, 2.0	21.0
482	Dissecting	1.0, 2.0	24.1
484	Gerald 1x2	-	17.8
486	Adson Bayonet	1.0, 2.0	16.0
488	Adson Bayonet	1.0, 2.0	19.1
490	Adson Bayonet Angled UP	1.0, 2.0	16.0
492	Adson Bayonet Angled UP	1.0, 2.0	19.1
494	Adson Bayonet Angled UP	1.0, 2.0	21.0
496	Remington Hobb	1.0, 2.0	24.1

ENGLISH MONOPOLAR MALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
496	Adson	1.0	12.0
498	McINDOE	1.0, 2.0	15.0
500	Gillies	-	15.0
502	Gillies	-	17.8
504	Gillies	-	20.0
506	Gillies	-	24.0
508	Dissecting	1.0, 2.0	15.0
510	Dissecting	1.0, 2.0	17.8
512	Dissecting	1.0, 2.0	21.0
514	Dissecting	1.0, 2.0	24.1
516	Gerald 1x2	-	17.8
518	Adson Bayonet	1.0, 2.0	16.0
520	Adson Bayonet	1.0, 2.0	19.1
522	Adson Bayonet Angled UP	1.0, 2.0	16.0
524	Adson Bayonet Angled UP	1.0, 2.0	19.1
526	Adson Bayonet Angled UP	1.0, 2.0	21.0
528	Remington Hobb	1.0, 2.0	24.1

MONOPOLAR WITHOUT CONNECTOR



Ref	Item Name	Tip Size	Length Size
530	Adson	1.0	12.0
532	McINDOE	1.0, 2.0	15.0
534	Gillies	-	15.0
536	Gillies	-	17.8
538	Gillies	-	20.0
540	Gillies	-	24.0
542	Dissecting	1.0, 2.0	15.0
544	Dissecting	1.0, 2.0	17.8
546	Dissecting	1.0, 2.0	21.0
548	Dissecting	1.0, 2.0	24.1
550	Gerald 1x2	-	17.8
552	Adson Bayonet	1.0, 2.0	16.0
554	Adson Bayonet	1.0, 2.0	19.1
556	Adson Bayonet Angled UP	1.0, 2.0	16.0
558	Adson Bayonet Angled UP	1.0, 2.0	19.1
560	Adson Bayonet Angled UP	1.0, 2.0	21.0
562	Remington Hobb	1.0, 2.0	24.1

Pro-Pas Electro disposable bipolar forceps have been designed and tested to meet the practical. The main purpose for offering single-use bipolar forceps is to avoid the risk of contamination while remaining cost effective and utilizing the functionalities of the higher end forceps. Hardened stainless steel provides an optimum spring tension and serves as the basis for a superior instruments.


Advantages of Disposable Bipolar Forceps

- 1- Lightweight and well balanced to prevent fatigue.*
- 2- Highly polished tips minimise sticking.*
- 3- Nylon coating resists abrasion and provides superb insulation.*
- 4- Precision tip alignment aids fine control.*
- 5- Cost effective in comparsion to re-usable.*
- 6- All Bipolar forceps available with cable, sterile packed.*
- 7- Compatible with all common electrosurgery unit.*







Ref	Item Name	Tip Size (mm)	Length Size (cm)
524	JEWELLER# 5	NEEDLE, 0.5	11.5
526	ADSON	0.5, 1.0	12
528	STRAIGHT	0.5, 1.0, 1.5	18
530	STRAIGHT ANGLED	0.5, 1.0, 1.5	18
532	BAYONET	0.5, 1.0, 1.5, 2.0	20
534	BAYONET ANGLED DOWN	0.5, 1.0, 1.5, 2.0	20




Ref	Item Name	Tip Size (mm)	Length Size (cm)
536	ADSON	0.5, 1.0	12
538	CUSHING	0.5, 1.0, 1.5	18
540	JANSEN BAYONET	0.5, 1.0, 1.5, 2.0	20
542	ADSON	0.5, 1.0	15
544	STRAIGHT	0.5, 1.0, 1.5, 2.0	20
546	STRAIGHT ANGLED	0.5, 1.0, 1.5, 2.0	20

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.






Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
548	BAYONET	0.5, 1.0, 1.5, 2.0	20
552	JEWELLER# 5	NEEDLE, 0.5	11.5
556	ADSON	0.5, 1.0	15
558	STRAIGHT	0.5, 1.0, 1.5, 2.0	20
560	STRAIGHT ANGLED	0.5, 1.0, 1.5, 2.0	20
562	BAYONET	0.5, 1.0, 1.5, 2.0	20


Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.


DISPOSABLE BIPOLAR FORCEPS



Ref	Item Name	Tip Size (mm)	Length Size (cm)
564	MCPHERSON	NEEDLE, 0.5	10.2
566	ADSON	0.5, 1.0	15
568	JEWELLER# 5	NEEDLE, 0.5	11.5
570	STRAIGHT	0.5, 1.0, 1.5, 2.0	20
574	BAYONET	0.5, 1.0, 1.5, 2.0	20
576	JEWELLER# 5 NON-STICK	NEEDLE, 0.5	11.5

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.





Ref	Item Name	Tip Size (mm)	Length Size (cm)
578	ADSON NON-STICK	0.5, 1.0	15
580	STRAIGHT NON-STICK	0.5, 1.0, 1.5, 2.0	20
582	BAYONET NON-STICK	0.5, 1.0, 1.5, 2.0	20

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.

The First Choice For Proven Performance, Pro-Pas Electro Provides The Broadest selection of fully autoclavable instruments for LEEP procedures. Our "rose coated" instruments will provide you with durability and performance unmatched in the industry. Models Available: Graves, Pederson, VU-More, 4-Way, Cusco, Collin, Weisman, Gelpi, Rigby, Lateral, Lateral Vaginal, Endocervical, Iris, Emmet Straight, two-prong, three-prong, etc.





Disposable Tube



Graves Vaginal Speculum

With Smoke Evacuation Tube

Small **381-600**

Medium **381-602**

Large **381-604**

Small **381-606**

Medium **381-608**

Large **381-610**



Disposable Tube



Pederson Vaginal Speculum

With Smoke Evacuation Tube

Small **381-612**

Medium **381-614**

Large **381-616**

Small **381-618**

Medium **381-620**

Large **381-622**



Disposable Tube



VU-Max Grave Vaginal Speculum

With Smoke Evacuation Tube

Medium **381-624**

Large **381-626**

Medium **381-628**

Large **381-630**



Disposable Tube



VU-Max Pederson Vaginal Speculum

With Smoke Evacuation Tube

Medium **381-632**

Large **381-634**

Medium **381-636**

Large **381-638**



Collin Vaginal Speculum

With Smoke Evacuation Tube

Small **381-640**

Medium **381-642**

Large **381-644**



CUSCO Vaginal Speculum

With Smoke Evacuation Tube

Small **381-646**

Medium **381-648**

Large **381-650**



Disposable Tube



VU-More Graves Vaginal Speculum
With Smoke Evacuation Tube
Medium **381-652**
Large **381-654**



4-Blades Speculum
With Smoke Evacuation Tube
Medium **381-660**
Large **381-662**



Disposable Tube



Weisman Vaginal Speculum
Right Side Open
With Smoke Evacuation Tube
Medium **381-664**
Large **381-666**



Disposable Tube



Weisman Vaginal Speculum
Left Side Open
With Smoke Evacuation Tube
Medium **381-672**
Large **381-674**

Medium **381-676**
Large **381-678**



Lateral Retractor
Wide Open
381-680



Rigby Retractor
Wide Open
381-682



Gelpi Retractor
Wide Open
381-684



Kogan Speculum
3mm **381-686**
4mm **381-688**
5mm **381-690**



Lateral Retractor
21cm 22mmx75mm
381-692



**Schroeder
Tanaculum 24cm**
381-694



**Schroeder
Tanaculum 24cm**
381-696



**Single Tooth
Tanaculum 24cm**
381-698



**Campion
Forceps 22cm**
381-700



**Spong
Forceps 22cm**
381-702

Dressing Forceps
381-704 21cm
381-706 25cm



Tissue Forceps
381-708 21cm
381-710 25cm



Iris Hook
381-712 25cm



Emmet Hook
381-714 25cm



Emmet Hook Straight
381-716 25cm



Emmet Hook Two Prong
381-718 25cm



Emmet Hook Three Prong
381-720 25cm



These E-Z Clean Laparoscopic Electrodes are made from very high quality raw material which ensures hassle work performance at its user end. These E-Z Clean Laparoscopic Electrodes are widely finds its applications in various medical sectors.

Salient Features of our E-Z Clean Laparoscopic Electrodes

- 1- Provide excellent dissection and coagulation for a variety of surgical applications.*
- 2- High-tech PTFE coating for reduced eschar buildup.*
- 3- Reintroduction during surgery is minimized.*
- 4- Lower power settings for safer delivery of electrosurgery.*
- 5- Tips are cost-effective because they Save Time and Money.*



Without Suction



L-Hook Electrode 36cm
381-800



J-Hook Electrode 36cm
381-802



Needle Electrode 36cm
381-804



Spectula Electrode 36cm
381-806



Ball Electrode 36cm
381-808



With Suction



L-Hook Electrode 36cm
381-810



J-Hook Electrode 36cm
381-812



Needle Electrode 36cm
381-814



Spectula Electrode 36cm
381-816



Ball Electrode 36cm
381-818



Electrode handle with blade electrode for single-use. Sterile packaging in individual pouches. Ready to use. This handle is made of high quality material and fits into the surgeon's hand comfortably due to its ergonomic shape. The quality meets highest expectations and fulfills all required international safety regulations.

Advantages of Reusable European Irrigation Bipolar Forceps

- 1-** *Single-use electrode handle including blade electrode and 3 m cable (also available: 5 m)*
- 2-** *Individually packaged, sterile.*
- 3-** *Push buttons for cutting and coagulation.*
- 4-** *Unit side: 3-pin connector.*
- 5-** *Tips are cost-effective because they Save Time and Money.*
- 6-** *Suitable also for use with FRBE, Martin or Berchtold units, by addition of an adopter especially designed for such purpose.*





Rocker Switch Pencil (Single Use)
3m Cable 381-854



Foot Control Pencil (Single Use)
3m Cable 381-856
5m Cable 381-858



Finger Switch Pencil (Reusable)
3m Silicon Cable 381-860



Electrode Handle (Reusable)
3m Silicon Cable 381-862



Finger Switch Pencil (Single Use)
3m Cable 381-850
5m Cable 381-852



Blade Electrode

Disposable: 381-900-S Reusable: 381-900-R



Blade Electrode Crv

Disposable: 381-902-S Reusable: 381-902-R



Needle Electrode

Disposable: 381-904-S Reusable: 381-904-R



Needle Electrode Crv

Disposable: 381-906-S Reusable: 381-906-R



Ball Electrode 2mm

Disposable: 381-908-S Reusable: 381-908-R



Ball Electrode Crv 2mm

Disposable: 381-910-S Reusable: 381-910-R



Ball Electrode 3mm

Disposable: 381-912-S Reusable: 381-912-R



Ball Electrode 3mm Crv

Disposable: 381-914-S Reusable: 381-914-R



Ball Electrode 4mm

Disposable: 381-916-S Reusable: 381-916-R



Ball Electrode 4mm Crv

Disposable: 381-918-S Reusable: 381-918-R



Ball Electrode 5mm

Disposable: 381-920-S Reusable: 381-920-R



Ball Electrode 5mm Crv

Disposable: 381-922-S Reusable: 381-922-R



Blade Electrode Shaft 13cm

Disposable: 381-924-S Reusable: 381-924-R



Needle Electrode Shaft 13cm

Disposable: 381-926-S Reusable: 381-926-R



Ball Electrode 3mm Shaft 13cm

Disposable: 381-928-S Reusable: 381-928-R



Ball Electrode 5mm Shaft 13cm

Disposable: 381-930-S Reusable: 381-930-R



Conization Electrode

Disposable: 381-932-S Reusable: 381-932-R



Circle Loop Electrode 4mm

Disposable: 381-934-S Reusable: 381-934-R



Circle Loop Electrode 6mm

Disposable: 381-936-S Reusable: 381-936-R



Circle Loop Electrode 8mm

Disposable: 381-938-S Reusable: 381-938-R



Diamond Electrode

Disposable: 381-940-S Reusable: 381-940-R



Fine Needle Electrode

Disposable: 381-942-S Reusable: 381-942-R



Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	381-944-S	381-946-S	381-948-S	381-950-S	381-952-S
Reusable	381-944-R	381-946-R	381-948-R	381-950-R	381-952-R



Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	381-954-S	381-956-S	381-958-S	381-960-S	381-962-S
Reusable	381-954-R	381-956-R	381-958-R	381-960-R	381-962-R



Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	381-964-S	381-966-S	381-968-S	381-970-S	381-972-S	381-974-S
Reusable	381-964-R	381-966-R	381-968-R	381-970-R	381-972-R	381-974-R



Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	381-976-S	381-978-S	381-980-S	381-982-S	381-984-S	381-986-S
Reusable	381-976-R	381-978-R	381-980-R	381-982-R	381-984-R	381-986-R



Blade Electrode

Disposable: 381-1000-S Reusable: 381-1000-R



Blade Electrode Crv

Disposable: 381-1002-S Reusable: 381-1002-R



Needle Electrode

Disposable: 381-1004-S Reusable: 381-1004-R



Needle Electrode Crv

Disposable: 381-1006-S Reusable: 381-1006-R



Ball Electrode 2mm

Disposable: 381-1008-S Reusable: 381-1008-R



Ball Electrode Crv 2mm

Disposable: 381-1010-S Reusable: 381-1010-R



Ball Electrode 3mm

Disposable: 381-1012-S Reusable: 381-1012-R



Ball Electrode 3mm Crv

Disposable: 381-1014-S Reusable: 381-1014-R



Ball Electrode 4mm

Disposable: 381-1016-S Reusable: 381-1016-R



Ball Electrode 4mm Crv

Disposable: 381-1018-S Reusable: 381-1018-R



Ball Electrode 5mm

Disposable: 381-1020-S Reusable: 381-1020-R



Ball Electrode 5mm Crv

Disposable: 381-1022-S Reusable: 381-1022-R



Blade Electrode Shaft 13cm

Disposable: 381-1024-S Reusable: 381-1024-R



Needle Electrode Shaft 13cm

Disposable: 381-1026-S Reusable: 381-1026-R



Ball Electrode Shaft 3mm 13cm

Disposable: 381-1028-S Reusable: 381-1028-R



Ball Electrode 4mm Shaft 13cm

Disposable: 381-1030-S Reusable: 381-1030-R



Ball Electrode 5mm 13cm

Disposable: 381-1032-S Reusable: 381-1032-R



Circle Loop Electrode 4mm 0.7cm

Disposable: 381-1034-S Reusable: 381-1034-R



Circle Loop Electrode 6mm 0.7cm

Disposable: 381-936-S Reusable: 381-936-R



Circle Loop Electrode 8mm 0.7cm

Disposable: 381-1038-S Reusable: 381-1038-R



Circle Loop Electrode 10mm 0.7cm

Disposable: 381-940-S Reusable: 381-940-R



Plate Electrode Small

Disposable: 381-1042-S Reusable: 381-1042-R



Plate Electrode Large

Disposable: 381-1044-S Reusable: 381-1044-R



Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	381-1046-S	381-1048-S	381-1050-S	381-1052-S	381-1054-S
Reusable	381-1046-R	381-1048-R	381-1050-R	381-1052-R	381-1054-R



Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	381-1056-S	381-1058-S	381-1060-S	381-1062-S	381-1064-S
Reusable	381-1056-R	381-1058-R	381-1060-R	381-1062-R	381-1064-R



Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	381-1066-S	381-1068-S	381-1070-S	381-1072-S	381-1074-S	381-1076-S
Reusable	381-1066-R	381-1068-R	381-1070-R	381-1072-R	381-1074-R	381-1076-R



Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	381-1078-S	381-1080-S	381-1082-S	381-1084-S	381-1086-S	381-1088-S
Reusable	381-1078-R	381-1080-R	381-1082-R	381-1084-R	381-1086-R	381-1088-R

Pro-Pas Electro Cables are available in disposable and reusable grade. We offer high quality silicon cables even for our single-use instruments thus making sure that all cables are able to withstand electrical current above the standard requirements. Furthermore, numerous tests including continuity and leak test ensures that our cables will work under the toughest requirements. Multiples sizes and plugs are available as per the need of the customer. Our cables are durable, flexible and able to withstand sterilization processes if needed. By following international standards we ensure that our cables are compatible with most major ESU's. Our cables are available in lengths 3 m, 4 m or 5 m.





Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **381-1100**
Single Use 3m Cable **381-1102**



Bipolar Cable

Instruments Connector: **European Flat Plug**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **381-1104**
Single Use 3m Cable **381-1106**



Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**
Generator Connector: **Pin Distance 29mm**
Reusable Silicon 3m Cable **381-1112**
Single Use 3m Cable **381-1114**



Bipolar Cable

Instruments Connector: **European Flat Plug**
Generator Connector: **Pin Distance 29mm**
Reusable Silicon 3m Cable **381-1108**
Single Use 3m Cable **381-1110**



Bipolar Cable

Instruments Connector: **Bipolar Artery Sealer**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **381-1116**



European Bipolar Cable

Instruments Connector: **European Flat Plug**

Generator Connector: **Martin Type**

Reusable Silicon 3m Cable **381-1118**



US Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**

Generator Connector: **Martin Type**

Reusable Silicon 3m Cable **381-1120**



European Bipolar Cable

Instruments Connector: **European Flat Plug**

Generator Connector: **Erbe Type**

Reusable Silicon 3m Cable **381-1122**



US Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**

Generator Connector: **Erbe Type**

Reusable Silicon 3m Cable **381-1124**



Bipolar Cable

Instruments Connector: **Ethicon**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **381-1126**



Monopolar Cable

Instruments Connector: **Screw Type**

Generator Connector: **Standard Bovie 8mm Plug**

Reusable Silicon 3m Cable **381-1128**



Monopolar Cable

Instruments Connector: **4mm Female Fitting**

Generator Connector: **Standard Bovie 8mm Plug**

Reusable Silicon 3m Cable **381-1130**



Monopolar Endoscopic Cable

Instruments Connector: **4mm Female Fitting**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **381-1132**



Diathermy Cable

Instruments Connector: **4mm Male Fitting**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **381-1134**



Diathermy Cable

Instruments Connector: **4.8mm Female Fitting**

Generator Connector: **4mm Banana Plug**

Reusable Silicon 3m Cable **381-1136**



Diathermy Cable

Instruments Connector: **4.8mm Female Fitting**

Generator Connector: **Stan. Bovie 8mm Plug**

Reusable Silicon 3m Cable **381-1138**

- 381-1150 Artery Sealer 6.5"
- 381-1152 Artery Sealer 7.5"
- 381-1154 Artery Sealer 8.5"
- 381-1156 Artery Sealer 9.5"
- 381-1158 Artery Sealer 10.5"



- 381-1160 Non-Stick Artery Sealer 6.5"
- 381-1162 Non-Stick Artery Sealer 7.5"
- 381-1164 Non-Stick Artery Sealer 8.5"
- 381-1166 Non-Stick Artery Sealer 9.5"
- 381-1168 Non-Stick Artery Sealer 10.5"



381-1170 BIPOLAR SCISSOR 18CM
381-1172 BIPOLAR SCISSOR 20CM



381-1174 BIPOLAR SCISSOR 18CM
381-1176 BIPOLAR SCISSOR 20CM





ELECTRO SURGICAL



CATALOGUE

Supplier of Quality Electro Surgical Instruments

Address: Model Town, Ugoki Road, Sialkot Pakistan.

Website: www.propasintl.com.pk

Email: marketing4propas@gmail.com

info@propasintl.com.pk

www.propasintl.com.pk